

# Washer Tag with Washer Disinfectant Cleaning Efficiency Test Soil




Washer Tags are designed to be used as a tag to identify if the tray is washed or not.

- ★ Ready to use (single-use product).
- ★ The test provides efficient and reliable results.
- ★ It is easy to read and interpret.
- ★ Suitable for archiving.
- ★ An important advantage is that information about which set is in it can be read without mixing the tray.
- ★ Patented Design.

**PATENTED**

**C-1064**



PRODUCT SIZE: 140 X 30 mm



BEFORE WASHING PROCESS



AFTER WASHING PROCESS

REF NO	DESCRIPTION	PACKAGING
C-1064	Washer Tag with Washer Disinfectant Cleaning Efficiency Test (Pack of 100)	50 Packs / CARTON

